



Biological medicine      Scientific experiment      Inkjet printing




Food machinery      Hydroponics experiment

## BPH Peristaltic Pump

Product manual

A/0 Version



- 1 -

### BPH Selection Guides

#### Parts introductions

##### Characteristics

- ◆ Small size, compact structure, weight only 110g
- ◆ Easy disassembly and maintenance, quick replacement of pump pipes
- ◆ Unique design, low noise



B:Baines P:Peristaltic pump H:Serial number

BPH — DK B20 A

1


2

3

1




2



| size: ID×OD( mm) | Pipe joint | size: ID×OD( mm) | Pipe joint |
|------------------|------------|------------------|------------|
| 1x3 (S10)        | 1/16       | 1x3 (B10)        | 1/16       |
| 2x4 (S20)        | 3/32       | 2x4 (B20)        | 3/32       |
| 3x5 (S30)        | 3/32       |                  |            |

\*For the corresponding flow rate, see page 3.  
Material of pipe joint: PP (polypropylene)

3



- 2 -

### BPH Selection Guide

#### Introduction of Pump Pipe and Corresponding Flow Rate

##### ■ Introductions of pump tube

| Code /Material of pump tube | brief introduction   | Life of pump tube |
|-----------------------------|--|-------------------|
| <b>B</b><br>PharMed®BPT     | Ozone and UV radiation resistance, heat resistance, ozone resistance, Acid and alkali resistance, anti-aging and anti-oxidation, long service life. Working temperature range from -59℃ ~ 135℃ | ≥1000h            |
| <b>S</b><br>Silicone tube   | Low cost, weak acid and alkali resistance, chemical resistance is worse and the service life is lower than the three pumping pipes above.  | ≥200h             |

※ The tubular test data are obtained by continuous operation of BPH 3 rotor pump head at 250 rpm/min at room temperature of 20 °C pure water without pressure until the cracking in the pump pipe appears. In general, the slower the speed of pump head, the longer the service life of pump pipe. The actual service life is affected by factors such as medium, working condition, voltage and so on. The test data can be used as reference.

##### ■ Flow rate

| Tubing              |                                 | S  | S  | S  | BPT | BPT | BPT |
|---------------------|---------------------------------|----|----|----|-----|-----|-----|
| Flow rate<br>ml/min | 24V<br>DC motor<br>(DK)<br>0.2A | 14 | 56 | 96 | 12  | 53  | 83  |
|                     | 12V<br>DC motor<br>(DL)<br>0.4A | 14 | 56 | 96 | 12  | 53  | 83  |
|                     | 6V<br>DC motor<br>(DN)<br>0.4A  | 14 | 56 | 96 | 12  | 53  | 83  |

Ideal working conditions: ambient temperature 0-40 °C relative humidity < 80%.

Ideal working conditions: ambient temperature 0-40 °C relative humidity < 80%.  
The total length of the pump pipe is 135 mm. The inlet and outlet of the pump pipe are respectively exposed to 30 mm of the pump head.


※ 1: The above data are measured by pure water without pressure at room temperature and standard atmospheric pressure at 20 °C. Actually, according to different medium, the outlet pressure is different. There will be some errors in the flow rate due to assembly tolerance, etc

- 3 -


### BPH Selection Guide

#### Installation mode


##### ■ Installation mode A



##### ■ Opening size A



##### ■ Installation mode B




- 4 -


### BPH Selection Guide

#### Sketch Map

##### ■ Installation diagram (A)




##### ■ Installation diagram (B)



##### ■ Assembly drawings

1.The motor shaft is inserted into the central hole of the pump head along the arrow direction.

2.Pump head is pushed in the direction of arrow and clamped on the motor board



- 5 -